


JM Disposable Dental Needles


- Indicator dot on hub provides easy identification of the lancet bevel position for less insertion pain at the injection site
- Lancet bevel point at cartridge end prevents blockage of anesthetic and ensures easy flow of anesthetic solution
- Universal plastic hub fits most syringes
- Color coding for easy identification
- Economical



Lancet bevel prevents rubber chips from blocking needle bores

One dot needle

Designed for both Luer-Lock and screw-on type syringes



Ordering Information

(100) sterile, disposable dental needles are included in each box.

25 G Long (30 mm)	20-25GL
27 G Short (21 mm)	20-27GS
27 G Long (30 mm)	20-27GL
30 G Short (21 mm)	20-30GS
30 G Long (25 mm)	20-30GL

