

Thinking ahead. Focused on life.

# New Tips for AdvErL Evo



Design enhancements for added safety & treatment of peri-implantitis

The AdvErL Evo Er:YAG laser is effective on both hard and soft tissue and can be used for a wide variety of procedures such as: implant recovery, treatment of peri-implantitis, soft tissue curettage, sulcular debridement and removal of subgingival calculus, removal of highly inflamed edematous tissue and granulation tissue from bony defects, and much more. This unit has an extensive tip product line for a multitude of clinical situations with an ideal shape for periodontal and implant treatment.





J. Morita has just announced three additional tips featuring brand new technology.

**R600T**



The R600T emits the laser from the sides of the tip, especially useful for debridement of peri-implantitis – this feature is completely unique in the market.

**CS600F**

The CS600F offers improved water spray for more effective cutting of enamel and dentin.

**PSM600T**

The PSM600T features a reinforced Titanium Nitride metal coating to reduce the risk of tip fracture and to make it easier to control inside the pocket.

## Ordering Information

R600T ..... 34-8924708  
 CS600F ..... 34-8924694  
 PSM600T ..... 34-8924716