



Tri Auto ZX2 Press Release

New product
May 2017

Tri Auto ZX2 Ingeniously versatile

Introducing the new Tri Auto ZX2, engineered with innovative and dynamic design form. As the successor to the original Tri Auto ZX cordless endodontic motor with integrated apex locator, Morita continues to offer the only endodontic system in the market combining both functions in one, cordless handpiece. With easy, intuitive operation and automatic functions, Tri Auto ZX2 delivers predictable and safe endodontic instrumentation.

Incorporating Morita's new OTR (Optimum Torque Reverse) and OGP (Optimum Glide Path) modes, Tri Auto ZX2 delivers, unique, compelling and reliable clinical application. Tri Auto ZX2: The complete endodontic preparation solution.

Exclusively available through Henry Schein Halas (Australia) and Henry Schein Dental (New Zealand).